Rhegnova





# Testing Sputum in Cystic Fibrosis Patients using Rheology: the Ageing Effect

Étienne Ghiringhelli,<sup>1</sup> Matthieu Robert de Saint Vincent,<sup>1</sup> Boubou Camara,<sup>2,3</sup> Sébastien Quétant,<sup>2,3</sup> Christophe Pison,<sup>2,3</sup> and Jérémy Patarin<sup>1</sup>

1 Rheonova, Grenoble, France, contact@rheonova.fr; 2 CHU Grenoble Alpes, Pôle Thorax et Vaisseaux, Clinique Universitaire de Pneumologie, Grenoble, France; 3 Univ. Grenoble-Alpes, LBFA-U1055, Grenoble, France

#### Context

**Sputum rheology** investigates mechanical properties of the mucus, *i.e.* elasticity and viscosity. Rheology is proposed as a prospective tool to **follow the status of patients with cystic fibrosis** (CF), and to **assess the individual efficiency of current and future treatments** during clinical practice (endpoints). However, this characterisation may involve sample ageing during their transient storage from collection to measurement. To better define rheological assessment protocols, the possible effect of **ageing** of sputum samples on their rheological properties is investigated.



**Sample Collection** 

#### **Sample Preparation**

Sputum is composed of highly viscoelastic mucus plugs included in a less viscoelastic matrix. This heterogeneity induces artefacts in rheology measurements.





### **Rheological Measurements**

**Rheology** is the study of the flow of matter that exhibits a combination of **elastic**, viscous and plastic behaviours by combining elasticity and fluid mechanics.



Results



## **Concluding Remarks**

**1** Mucus globally thickens when stored at 4°C. This thickening effect mainly occurs within the first 4–6 h after sputum collection. However, no significant gelation or fluidisation is observed, suggesting that this thickening is not related to evaporation.

G' and G'' both increase in **similar proportions** (x2 within 8-10 h).

 $\geq$  tan  $\delta$  remains constant.

2 Since the sputum viscoelasticity evolves with storage time, care must be taken to minimise it and/or to define unified protocols if rheology is to be used as a clinical marker.

#### Perspectives

Measurements performed at 4°C mimic the standard sample storage conditions. Other relevant storage conditions (*e.g.*, at room temperature (20°C), within controlled humidity, after freezing/thawing cycles) would call for complementary investigations to understand the possible changes in the mucus microstructure.